

INPP5K

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 51kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human INPP5K

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PPS; SKIP;

Catalog #: A7807

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51763

Isotype: IgG

Swiss Prot: Q9BT40

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein with 5-phosphatase activity toward polyphosphate inositol. The protein localizes to the cytosol in regions lacking actin stress fibers. It is thought that this protein may negatively regulate the actin cytoskeleton. Alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).