

SLC22A11

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 60kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SLC22A11

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OAT4; hOAT4;

Catalog #: A7816

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55867

Isotype: IgG

Swiss Prot: Q9NSA0

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is involved in the sodium-independent transport and excretion of organic anions, some of which are potentially toxic. The encoded protein is an integral membrane protein and is found mainly in the kidney and in the placenta, where it may act to prevent potentially harmful organic anions from reaching the fetus. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).