

ALG1

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 53kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ALG1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HMT1; MT-1; CDG1K; HMAT1; HMT-1; Mat-1; hMat-1;

Background:

The enzyme encoded by this gene catalyzes the first mannosylation step in the biosynthesis of lipid-linked oligosaccharides. This gene is mutated in congenital disorder of glycosylation type Ik.

To place an order, please [Click HERE](#).

Catalog #: A7818

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 56052

Isotype: IgG

Swiss Prot: Q9BT22

Purity: Affinity purification

For research use only.