

REEP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 22kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human REEP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HMN5B; SPG31; C2orf23;

Catalog #: A7832

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 65055

Isotype: IgG

Swiss Prot: Q9H902

Purity: Affinity purification

For research use only.

Background:

This gene encodes a mitochondrial protein that functions to enhance the cell surface expression of odorant receptors. Mutations in this gene cause spastic paraplegia autosomal dominant type 31, a neurodegenerative disorder. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).