

LRG1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:38kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human LRG1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LRG; HMFT1766;

Catalog #:A7850

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:116844

Isotype:IgG

Swiss Prot:P02750

Purity:Affinity purification

For research use only.

Background:

The leucine-rich repeat (LRR) family of proteins, including LRG1, have been shown to be involved in protein-protein interaction, signal transduction, and cell adhesion and development. LRG1 is expressed during granulocyte differentiation (O'Donnell et al., 2002)

To place an order, please [Click HERE](#).