

CALHM1

Reactivity:Human Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:38kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CALHM1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FAM26C;

Catalog #:A7858

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:255022

Isotype:IgG

Swiss Prot:Q8IU99

Purity:Affinity purification

For research use only.

Background:

This gene encodes a calcium channel that plays a role in processing of amyloid-beta precursor protein. A polymorphism at this locus has been reported to be associated with susceptibility to late-onset Alzheimer's disease in some populations, but the pathogenicity of this polymorphism is unclear.

To place an order, please [Click HERE](#).