

## ELMOD2

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 35kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human ELMOD2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

b

**Synonym:**

9830169G11Rik;

**Catalog #:** A7859

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 255520

**Isotype:** IgG

**Swiss Prot:** Q8IZ81

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes one of six engulfment and motility (ELMO) domain-containing proteins. This gene is thought to play a role in antiviral responses. Mutations in this gene may be involved in the cause of familial idiopathic pulmonary fibrosis.

*To place an order, please [Click HERE](#).*