

TSSK4

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 37kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TSSK4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TSK4; TSK-4; TSSK5; STK22E; C14orf20;

Catalog #: A7861

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 283629

Isotype: IgG

Swiss Prot: Q6SA08

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the testis-specific serine/threonine kinase family. The encoded protein is thought to be involved in spermatogenesis via stimulation of the CREB/CRE responsive pathway through phosphorylation of the cAMP responsive element binding protein transcription factor. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).