

ADCY3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:200

Calculated MW: 129kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ADCY3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AC3;

Catalog #: A7870

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 109

Isotype: IgG

Swiss Prot: O60266

Purity: Affinity purification

For research use only.

Background:

This gene encodes adenylyl cyclase 3 which is a membrane-associated enzyme and catalyzes the formation of the secondary messenger cyclic adenosine monophosphate (cAMP). This protein appears to be widely expressed in various human tissues and may be involved in a number of physiological and pathophysiological metabolic processes.

To place an order, please [Click HERE](#).