

## CHD1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:100 - 1:200

**Calculated MW:**197kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human CHD1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

f

**Catalog #:**A7883

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**1105

**Isotype:**IgG

**Swiss Prot:**O14646

**Purity:**Affinity purification

For research use only.

**Background:**

The CHD family of proteins is characterized by the presence of chromo (chromatin organization modifier) domains and SNF2-related helicase/ATPase domains. CHD genes alter gene expression possibly by modification of chromatin structure thus altering access of the transcriptional apparatus to its chromosomal DNA template.

*To place an order, please [Click HERE](#).*