

MEF2A

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50- 1:200

Calculated MW:55kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human MEF2A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

mef2; ADCAD1; RSRFC4; RSRFC9

Catalog #:A7911

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4205

Isotype:IgG

Swiss Prot:Q02078

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a DNA-binding transcription factor that activates many muscle-specific, growth factor-induced, and stress-induced genes. The encoded protein can act as a homodimer or as a heterodimer and is involved in several cellular processes, including muscle development, neuronal differentiation, cell growth control, and apoptosis. Defects in this gene could be a cause of autosomal dominant coronary artery disease 1 with myocardial infarction (ADCAD1). Several transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).