

STX1A

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 33kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human STX1A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

STX1; HPC-1; P35-1; SYN1A;

Catalog #: A7931

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6804

Isotype: IgG

Swiss Prot: Q16623

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the syntaxin superfamily. Syntaxins are nervous system-specific proteins implicated in the docking of synaptic vesicles with the presynaptic plasma membrane.

Syntaxins possess a single C-terminal transmembrane domain, a SNARE

To place an order, please [Click HERE](#).