

WNT9A

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 40kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human WNT9A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

WNT14;

Catalog #: A7939

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7483

Isotype: IgG

Swiss Prot: O14904

Purity: Affinity purification

For research use only.

Background:

The WNT gene family consists of structurally related genes that encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It is expressed in gastric cancer cell lines. The protein encoded by this gene shows 75% amino acid identity to chicken Wnt14, which has been shown to play a central role in initiating synovial joint formation in the chick limb. This gene is clustered with another family member, WNT3A, in the chromosome 1q42 region.

To place an order, please [Click HERE](#).