

PIK3R3

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:54kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PIK3R3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

opqr

Synonym:

p55; p55PIK; p55-GAMMA;

To place an order, please [Click HERE](#).

Catalog #:A7942

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8503

Isotype:IgG

Swiss Prot:Q92569

Purity:Affinity purification

For research use only.