

DIRAS3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 26kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human DIRAS3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ARHI; NOEY2;

Catalog #: A7950

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9077

Isotype: IgG

Swiss Prot: O95661

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the ras superfamily, and is expressed in normal ovarian and breast epithelial cells, but not in ovarian and breast cancers. It is an imprinted gene, with monoallelic expression of the paternal allele, which is associated with growth suppression. Thus, this gene appears to be a putative tumor suppressor gene whose function is abrogated in ovarian and breast cancers.

To place an order, please [Click HERE](#).