

PRDM7

Reactivity:Human Mouse Rat

Tested applications:WB IF

Recommended Dilution:WB 1:500 - 1:2000 IF 1:50 - 1:100

Calculated MW:56kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PRDM7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PFM4; ZNF910;

Background:

This gene encodes a member of a family of proteins that may have roles in transcription and other nuclear processes. The encoded protein contains a KRAB (Kruppel-associated box) domain -A box and a SET (Su(var)3-9, Enhancer-of-zeste, Trithorax) domain and may function as a histone methyltransferase.

To place an order, please [Click HERE](#).

Catalog #:A7961

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:11105

Isotype:IgG

Swiss Prot:Q9NQW5

Purity:Affinity purification

For research use only.