

## FTSJ1

---

**Reactivity:**Human

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**36kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human FTSJ1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

JM23; MRX9; SPB1; CDLIV; MRX44; TRMT7;

**Background:**

This gene encodes a member of the methyltransferase superfamily. The encoded protein localizes to the nucleolus, binds to S-adenosylmethionine, and may be involved in the processing and modification of ribosomal RNA. Mutations in this gene are associated with mental retardation. Alternative splicing results in multiple transcript variants.

*To place an order, please [Click HERE](#).*

**Catalog #:**A7967

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**24140

**Isotype:**IgG

**Swiss Prot:**Q9UET6

**Purity:**Affinity purification

For research use only.