

NPTN

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:44kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human NPTN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GP55; GP65; SDR1; np55; np65; SDFR1;

Background:

This gene encodes a type I transmembrane protein belonging to the Ig superfamily. The protein is believed to be involved in cell-cell interactions or cell-substrate interactions. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).

Catalog #:A7972

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:27020

Isotype:IgG

Swiss Prot:Q9Y639

Purity:Affinity purification

For research use only.