

YARS2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 53kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human YARS2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

h

Synonym:

TYRRS; CGI-04; MLASA2; MT-TYRRS;

Catalog #: A7976

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51067

Isotype: IgG

Swiss Prot: Q9Y2Z4

Purity: Affinity purification

For research use only.

Background:

This gene encodes a mitochondrial protein that catalyzes the attachment of tyrosine to tRNA(Tyr). Mutations in this gene are associated with myopathy with lactic acidosis and sideroblastic anemia type 2 (MLASA2).

To place an order, please [Click HERE](#).