

## IL20RB

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**35kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human IL20RB

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

DIRS1; FNDC6; IL-20R2;

**Background:**

IL20RB and IL20RA (MIM 605620) form a heterodimeric receptor for interleukin-20 (IL20; MIM 605619) (Blumberg et al., 2001)

*To place an order, please [Click HERE](#).*

**Catalog #:**A7980

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**53833

**Isotype:**IgG

**Swiss Prot:**Q6UXL0

**Purity:**Affinity purification

For research use only.