

RNF125

Reactivity:Human

Tested applications:WB IF

Recommended Dilution:WB 1:500 - 1:2000 IF 1:50 - 1:100

Calculated MW:26kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RNF125

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TNORS; TRAC1; TRAC-1;

Background:

This gene encodes a novel E3 ubiquitin ligase that contains a RING finger domain in the N-terminus and three zinc-binding and one ubiquitin-interacting motif in the C-terminus. As a result of myristoylation, this protein associates with membranes and is primarily localized to intracellular membrane systems. The encoded protein may function as a positive regulator in the T-cell receptor signaling pathway.

To place an order, please [Click HERE](#).

Catalog #:A7982

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:54941

Isotype:IgG

Swiss Prot:Q96EQ8

Purity:Affinity purification

For research use only.