

OPA3

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 20kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human OPA3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MGA3;

Catalog #: A7997

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 80207

Isotype: IgG

Swiss Prot: Q9H6K4

Purity: Affinity purification

For research use only.

Background:

The mouse ortholog of this protein co-purifies with the mitochondrial inner membrane. Mutations in this gene have been shown to result in 3-methylglutaconic aciduria type III and autosomal dominant optic atrophy and cataract. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).