

## DOCK8

**Reactivity:** Human Mouse Rat

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 239kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human DOCK8

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

MRD2; ZIR8; HEL-205;

**Catalog #:** A7998

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 81704

**Isotype:** IgG

**Swiss Prot:** Q8NF50

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a member of the DOCK180 family of guanine nucleotide exchange factors. Guanine nucleotide exchange factors interact with Rho GTPases and are components of intracellular signaling networks. Mutations in this gene result in the autosomal recessive form of the hyper-IgE syndrome. Alternatively spliced transcript variants encoding different isoforms have been described.

*To place an order, please [Click HERE](#).*