

USP26

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:104kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human USP26

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A7999

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:83844

Isotype:IgG

Swiss Prot:Q9BXU7

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the ubiquitin-specific processing (UBP) family of proteases and is a deubiquitinating enzyme (DUB) with His and Cys domains. It is specifically expressed in testis tissue. Mutations in this gene have been associated with Sertoli cell-only syndrome and male infertility.

To place an order, please [Click HERE](#).