

EGLN3

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:27kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human EGLN3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bf

Synonym:

PHD3; HIFPH3; HIFP4H3;

To place an order, please [Click HERE](#).

Catalog #:A8001

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:112399

Isotype:IgG

Swiss Prot:Q9H6Z9

Purity:Affinity purification

For research use only.