

## EFHC1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**74kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human EFHC1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

EJM1; dJ304B14.2;

**Catalog #:**A8002

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**114327

**Isotype:**IgG

**Swiss Prot:**Q5JVL4

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes an EF-hand-containing calcium binding protein. The encoded protein likely plays a role in calcium homeostasis. Mutations in this gene have been associated with susceptibility to juvenile myoclonic epilepsy and juvenile absence epilepsy. Alternatively spliced transcript variants have been described.

*To place an order, please [Click HERE](#).*