

CEP120

Reactivity:Human Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:113kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CEP120

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SRTD13; CCDC100;

Catalog #:A8004

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:153241

Isotype:IgG

Swiss Prot:Q8N960

Purity:Affinity purification

For research use only.

Background:

This gene encodes a protein that functions in the microtubule-dependent coupling of the nucleus and the centrosome. A similar protein in mouse plays a role in both interkinetic nuclear migration, which is a characteristic pattern of nuclear movement in neural progenitors, and in neural progenitor self-renewal. Mutations in this gene are predicted to result in neurogenic defects. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).