

NDUFS3

Reactivity: Human Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 30kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NDUFS3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bi

Synonym:

CI-30;

Catalog #: A8013

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4722

Isotype: IgG

Swiss Prot: O75489

Purity: Affinity purification

For research use only.

Background:

This gene encodes one of the iron-sulfur protein (IP) components of mitochondrial NADH:ubiquinone oxidoreductase (complex I). Mutations in this gene are associated with Leigh syndrome resulting from mitochondrial complex I deficiency.

To place an order, please [Click HERE](#).