

ITGA6

Reactivity: Human

Tested applications: IHC

Recommended Dilution: IHC 1:50 - 1:200

Calculated MW: 127kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human ITGA6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD49f; VLA-6; ITGA6B;

Catalog #: A8025

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3655

Isotype: IgG

Swiss Prot: P23229

Purity: Affinity purification

For research use only.

Background:

The ITGA6 protein product is the integrin alpha chain alpha 6. Integrins are integral cell-surface proteins composed of an alpha chain and a beta chain. A given chain may combine with multiple partners resulting in different integrins. For example, alpha 6 may combine with beta 4 in the integrin referred to as TSP180, or with beta 1 in the integrin VLA-6. Integrins are known to participate in cell adhesion as well as cell-surface mediated signalling. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).