

SSR4

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:19kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SSR4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SSR delta, SSR4, TRAP delta, TRAPD;

Catalog #:A8037

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6748

Isotype:IgG

Swiss Prot:P51571

Purity:Affinity purification

For research use only.

Background:

Signal sequence receptor subunit delta (SSR4, or TRAPD) gene encodes the translocon-associated proteins subunit delta, which is about 19-kDa and is located at endoplasmic reticulum (ER) membrane, the translocon-associated proteins complex functions to bind calcium to the ER membrane and thereby regulate the retention of ER resident proteins. TRAP-delta can form heterotetramer with other TRAP homolog TRAP-alpha, TRAP-beta and TRAP-gamma.

To place an order, please [Click HERE](#).