

BLVRA

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 33kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human BLVRA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BVR; BLVR; BVRA;

Catalog #: A8046

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 644

Isotype: IgG

Swiss Prot: P53004

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the biliverdin reductase family, members of which catalyze the conversion of biliverdin to bilirubin in the presence of NADPH or NADH. Mutations in this gene are associated with hyperbiliverdinemia. Alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).