

MTCH1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:42kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human MTCH1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PSAP; PIG60; CGI-64; SLC25A49;

Background:

This gene encodes a member of the mitochondrial carrier family. The encoded protein is localized to the mitochondrion inner membrane and induces apoptosis independent of the proapoptotic proteins Bax and Bak. Pseudogenes on chromosomes 6 and 11 have been identified for this gene. Alternatively spliced transcript variants encoding multiple isoforms have been observed.

To place an order, please [Click HERE](#).

Catalog #:A8063

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23787

Isotype:IgG

Swiss Prot:Q9NZJ7

Purity:Affinity purification

For research use only.