

TNFAIP8

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:23kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TNFAIP8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

NDED; GG2-1; SCCS2; SCC-S2; MDC-3.13;

To place an order, please [Click HERE](#).

Catalog #:A8064

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:25816

Isotype:IgG

Swiss Prot:O95379

Purity:Affinity purification

For research use only.