

PHF19

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:66kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PHF19

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PCL3; MTF2L1; TDRD19B;

To place an order, please [Click HERE](#).

Catalog #:A8065

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:26147

Isotype:IgG

Swiss Prot:Q5T6S3

Purity:Affinity purification

For research use only.