

ITGBL1

Reactivity:Human

Tested applications:WB IHC FC

Recommended Dilution:WB 1:200 - 1:1000 IHC 1:50 - 1:100 FC 1:10 - 1:50

Calculated MW:54kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human ITGBL1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OSCP; TIED;

Catalog #:A8084

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9358

Isotype:IgG

Swiss Prot:O95965

Purity:Affinity purification

For research use only.

Background:

This gene encodes a beta integrin-related protein that is a member of the EGF-like protein family.

The encoded protein contains integrin-like cysteine-rich repeats. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).