

## MonoMethyl-NONO-R357/R365

**Reactivity:** Human

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 54kDa

**Observed MW:** Refer to figures

**Immunogen:**

A synthetic peptide of human NONO

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

j

**Synonym:**

P54; NMT55; NRB54; P54NRB; PPP1R114;

**Catalog #:** A8087

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 4841

**Isotype:** IgG

**Swiss Prot:** Q15233

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes an RNA-binding protein which plays various roles in the nucleus, including transcriptional regulation and RNA splicing. A rearrangement between this gene and the transcription factor E3 gene has been observed in papillary renal cell carcinoma. Alternatively spliced transcript variants have been described. Pseudogenes exist on Chromosomes 2 and 16.

*To place an order, please [Click HERE](#).*