

RAB3IP

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:53kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human RAB3IP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RABIN3; RABIN8

Background:

RAB3IP (RAB3A interacting protein) is a protein-coding gene. GO annotations related to this gene include Rab guanyl-nucleotide exchange factor activity. An important paralog of this gene is RAB31L1.

To place an order, please [Click HERE](#).

Catalog #:A8094

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:117177

Isotype:IgG

Swiss Prot:Q96QF0

Purity:Affinity purification

For research use only.