

Gm13125

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW:55kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human Gm13125

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

EG627009; OTTMUSG00000010547;

To place an order, please [Click HERE](#).

Catalog #:A8095

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:645425

Isotype:IgG

Swiss Prot:B1ARV6

Purity:Affinity purification

For research use only.