

IL33

Reactivity: Human Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 31kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IL33

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DVS27; IL1F11; NF-HEV; NFEHEV; C9orf26;

Catalog #: A8096

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 90865

Isotype: IgG

Swiss Prot: O95760

Purity: Affinity purification

For research use only.

Background:

IL33 (MIM 608678) is a member of the IL1 (see MIM 147760) family that potently drives production of T helper-2 (Th2)-associated cytokines (e.g., IL4; MIM 147780). IL33 is a ligand for IL33R (IL1RL1; MIM 601203), an IL1 family receptor that is selectively expressed on Th2 cells and mast cells (summary by Yagami et al., 2010)

To place an order, please [Click HERE](#).