

CDH6

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 88kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CDH6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CAD6; KCAD;

Catalog #: A8109

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1004

Isotype: IgG

Swiss Prot: P55285

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the cadherin superfamily. Cadherins are membrane glycoproteins that mediate homophilic cell-cell adhesion and play critical roles in cell differentiation and morphogenesis. The encoded protein is a type II cadherin and may play a role in kidney development as well as endometrium and placenta formation. Decreased expression of this gene may be associated with tumor growth and metastasis.

To place an order, please [Click HERE](#).