

CEACAM7

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:29kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CEACAM7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CGM2;

Catalog #:A8112

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1087

Isotype:IgG

Swiss Prot:Q14002

Purity:Affinity purification

For research use only.

Background:

This gene encodes a cell surface glycoprotein and member of the carcinoembryonic antigen (CEA) family of proteins. Expression of this gene may be downregulated in colon and rectal cancer. Additionally, lower expression levels of this gene may be predictive of rectal cancer recurrence. This gene is present in a CEA family gene cluster on chromosome 19. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).