

CSTF2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 61kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CSTF2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CstF-64;

Catalog #: A8116

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1478

Isotype: IgG

Swiss Prot: P33240

Purity: Affinity purification

For research use only.

Background:

This gene encodes a nuclear protein with an RRM (RNA recognition motif) domain. The protein is a member of the cleavage stimulation factor (CSTF) complex that is involved in the 3' end cleavage and polyadenylation of pre-mRNAs. Specifically, this protein binds GU-rich elements within the 3'-untranslated region of mRNAs.

To place an order, please [Click HERE](#).