

ERCC4

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:104kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ERCC4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

XPF; RAD1; FANCD1; ERCC11;

Catalog #:A8119

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2072

Isotype:IgG

Swiss Prot:Q92889

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene forms a complex with ERCC1 and is involved in the 5' incision made during nucleotide excision repair. This complex is a structure specific DNA repair endonuclease that interacts with EME1. Defects in this gene are a cause of xeroderma pigmentosum complementation group F (XP-F), or xeroderma pigmentosum VI (XP6).

To place an order, please [Click HERE](#).