

NEU2

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:42kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human NEU2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SIAL2;

Catalog #:A8137

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4759

Isotype:IgG

Swiss Prot:Q9Y3R4

Purity:Affinity purification

For research use only.

Background:

This gene belongs to a family of glycohydrolytic enzymes which remove sialic acid residues from glycoproteins and glycolipids. Expression studies in COS7 cells confirmed that this gene encodes a functional sialidase. Its cytosolic localization was demonstrated by cell fractionation experiments.

To place an order, please [Click HERE](#).