

SRM

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:34kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SRM

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PAPT; SPS1; SPDSY; SRML1;

Catalog #:A8151

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6723

Isotype:IgG

Swiss Prot:P19623

Purity:Affinity purification

For research use only.

Background:

The polyamines putrescine, spermine, and spermidine are ubiquitous polycationic mediators of cell growth and differentiation. Spermidine synthase is one of four enzymes in the polyamine-biosynthetic pathway and carries out the final step of spermidine biosynthesis. This enzyme catalyzes the conversion of putrescine to spermidine using decarboxylated S-adenosylmethionine as the cofactor.

To place an order, please [Click HERE](#).