

SLC25A11

Reactivity:Human

Tested applications:WB IP

Recommended Dilution:WB 1:500 - 1:2000 IP 1:20 - 1:100

Calculated MW:34kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SLC25A11

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OGC; SLC20A4;

Background:

The oxoglutarate/malate carrier transports 2-oxoglutarate across the inner membranes of mitochondria in an electroneutral exchange for malate or other dicarboxylic acids (summary by Iacobazzi et al., 1992)

To place an order, please [Click HERE](#).

Catalog #:A8163

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8402

Isotype:IgG

Swiss Prot:Q02978

Purity:Affinity purification

For research use only.