

TAF1B

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 69kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TAF1B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SL1; RAF1B; RAFI63; TAFI63; MGC:9349;

Catalog #: A8171

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9014

Isotype: IgG

Swiss Prot: Q53T94

Purity: Affinity purification

For research use only.

Background:

Initiation of transcription by RNA polymerase I requires the formation of a complex composed of the TATA-binding protein (TBP) and three TBP-associated factors (TAFs) specific for RNA polymerase I. This complex, known as SL1, binds to the core promoter of ribosomal RNA genes to position the polymerase properly and acts as a channel for regulatory signals. This gene encodes one of the SL1-specific TAFs.

To place an order, please [Click HERE](#).