

CLDN8

Reactivity: Human Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 25kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CLDN8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A8174

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9073

Isotype: IgG

Swiss Prot: P56748

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the claudin family. Claudins are integral membrane proteins and components of tight junction strands. Tight junction strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epithelial or endothelial cell sheets, and also play critical roles in maintaining cell polarity and signal transductions. This protein plays important roles in the paracellular cation barrier of the distal renal tubule, and in the paracellular barrier to prevent sodium back-leakage in distal colon. Differential expression of this gene has been observed in colorectal carcinoma and renal cell tumors, and along with claudin-7, is an immunohistochemical marker for the differential diagnosis of chromophobe renal cell carcinoma and renal oncocytoma.

To place an order, please [Click HERE](#).