

METTL7A

Reactivity: Human Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 28kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human METTL7A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AAM-B;

To place an order, please [Click HERE](#).

Catalog #: A8201

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 25840

Isotype: IgG

Swiss Prot: Q9H8H3

Purity: Affinity purification

For research use only.