

IL36A

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:18kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human IL36A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FIL1; FIL1E; IL1F6; IL-1F6; IL1(EPSILON); FIL1(EPSILON);

To place an order, please [Click HERE](#).

Catalog #:A8207

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:27179

Isotype:IgG

Swiss Prot:Q9UHA7

Purity:Affinity purification

For research use only.